

## **PM Resting Respirometry Checklist**

- ☐ Gas tank volumes checked (should be >500 psi)
- ☐ System rinsed with DI water
- ☐ Chambers rinsed with hot water
- ☐ Chiller filled with packed ice/water and hoses not kinked
- ☐ Confirm MO2s don't decrease during measurement period with no fish (microbial respiration)
- ☐ Probes in correct lines

| Probe Label Color | Witrox Channel/Chamber # |
|-------------------|--------------------------|
| White             | Channel 1/ Chamber 1     |
| Red               | Channel 2/ Chamber 2     |
| Green             | Channel 3/ Chamber 3     |
| Blue              | Channel 4/ Chamber 4     |

- ☐ Water bath recirculation pump plugged in
  - ☐ N2 and CO2 gas flow adjusted and monitored
  - ☐ N2 gas tank turned off for night if running Control or 7.5 pH fish
  - ☐ pH/O2 treatment levels confirmed with Hach probe
  - ☐ Water temperature at 12° C
  - ☐ Barometric pressure updated
  - ☐ Flush/Wait/Measure periods set to 300/60/300
  - ☐ "Recirc pump always on" box **UNCHECKED**
  - ☐ Ramp setpoint inputs ready for Pcrit and unchecked
  - ☐ Air bubbles removed
  - ☐ Flushing outflow rate visually checked and equal between chambers
  - ☐ Hose clamps checked (don't adjust clamps unless something looks off)
  - ☐ Double check fish IDs and weight inputs in AutoResp
  - ☐ Chamber and tube volumes correct (128/46) and in mL
  - ☐ Double check for air bubbles after loading fish
  - ☐ Chambers closed tightly
  - ☐ MO2 levels checked (should be 80-200 mgO2/kg/hr)
  - ☐ Chambers reaching correct and similar DO levels after flushing (within 1% saturation)
  - ☐ AutoResp set to Intermittent Flow
  - ☐ Gas tank psi re-checked
  - ☐ Lights off in respirometry area
- computer monitor off

## **AM Pcrit Respirometry Checklist**

**\*\* In case of an emergency (if the DO levels drop too far or overshoot), press “skip phase” to skip a measurement phase and return to flushing without stopping the experiment.**

- ☐ Check MO<sub>2</sub>, DO, pH, and temp levels in AutoResp
- ☐ Chiller refilled with ice
- ☐ Barometric pressure updated
- ☐ Water bath recirculation pump plugged in
- ☐ Do not quit AutoResp program without first choosing “End Experiment”

### **Control or 7.5 pH treatment:**

- ☐ New experiment started and run for 3 loops at ambient
- ☐ Carefully lower DO to 70% while adjusting N<sub>2</sub> flow rate during 4<sup>th</sup> flush period

### **4.0 DO or X treatment:**

- ☐ New experiment started and DO raised to 70% with air stone
- ☐ Acclimate 1 hour on intermittent flow
- ☐ Start new experiment

### **Ramping Setpoints:**

- ☐ Chamber DO's at 70- if lower, raise water bath setpoint to ~72
- ☐ DO levels not overshooting by more than 2%
- ☐ Ramp setpoint to 40 by -10% every 3 loops (takes about 1 hour 10 mins)
- ☐ Chambers rising to correct and similar DO during wait periods (within 1% saturation)
- ☐ O<sub>2</sub> interval changed to -5 and minimum to 15% during one of the three 40% loops (takes ~2 hr 45 mins to reach 15% loops)
- ☐ Chiller ice refreshed
- ☐ System monitored during 15% loops – make sure chamber DOs reading above 8% during wait and measurement phase

### **Clean Up**

- ☐ After last loop at 15% take screenshot of MO<sub>2</sub> vs time and MO<sub>2</sub> vs DO graphs and save to Pictures -> OA/DO Respirometry folder
- ☐ Click “End Experiment”
- ☐ Remove fish as soon as possible, checking each tag to confirm chamber locations
- ☐ Trial data backed up to external hard drive
- ☐ Respirometry cooler drained, rinsed with DI
- ☐ Cooler, tubing and chambers dried thoroughly if not reloading
- ☐ Chiller ice drained